CUSTOMER NOTICE

Due to the properties of reagents used in IHC, allowances may be required for the environmental conditions found at each site. In particular, the boiling point of Epitope Retrieval solutions varies with air pressure as experienced at different altitudes. This means that sites at or above 250 mASL (820 ftASL) may experience a lower Epitope Retrieval solution boiling point which could lead to an area of suboptimal chromogenic staining.

Each site at the altitude ranges listed below should therefore consider their altitude when setting Epitope Retrieval (ER) temperatures as part of optimization and validation.

The table below provides suggested temperature settings for different altitude ranges for Heat Induced Epitope Retrieval (HIER) protocols. If your laboratory consistently achieves high quality stains using alternate temperature settings to those proposed in [Table 1-1](#), there is no need to adopt the suggested temperatures.

BOND-PRIME is not suitable for altitudes above 2250 meters Above Sea Level (mASL) (7381 feet Above Sea Level (ftASL)).

Table 1-1: Recommended HIER protocol temperature by altitude range

Altitude range (mASL)	Altitude range (ftASL)	Temperature in software (°C)
0 – 249	0 – 819	104 (default)
250 – 449	820 – 1475	102
450 – 799	1476 – 2624	101
800 – 1149	2625 – 3772	100
1150 – 1499	3773 – 4920	99
1500 – 1849	4921 – 6069	98
1850 – 2250	6070 – 7381	96

Contact your Leica Biosystems representative for a step-by-step guide on modifying the HIER temperature for BOND-PRIME. You can also refer to Chapter 7 of the *BOND 7 User Manual* for guidance on how to edit protocols.

Note that these are suggested temperatures only and each laboratory is responsible for validating assays in their environment as staining can be impacted by many factors. If you experience staining issues, contact your Leica Biosystems representative before making the changes to the predefined protocols.

The maximum altitude figure in this communication replaces the figure in the *BOND-PRIME User Manual (91.7500.501)*. The recommended temperature settings are in addition to the instructions in the *BOND 7 User Manual (49.7556.501)*.

For In Vitro Diagnostic Use

Leica and the Leica logo are registered trademarks of Leica Microsystems IR GmbH and used under license. BOND and BOND-PRIME are trademarks of Leica Biosystems Melbourne Pty Ltd ACN 008 582 401. Other trademarks are the property of their owners.

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Copyright © 2025 Leica Biosystems Melbourne Pty Ltd
© 2025 Leica Biosystems Melbourne Pty Ltd